

UBN2 Rabbit pAb

Catalog NO.: A10516 **1 Publications**

Basic Information

Observed MW

190kDa

Calculated MW

130kDa/131kDa/141kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Predicted to be involved in DNA replication-independent chromatin assembly. Predicted to be located in nucleoplasm. Predicted to be active in nucleus. Is expressed in central nervous system and retina. Orthologous to human UBN2 (ubinnuclein 2).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
320538	

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

mKIAA2030; 2900060J04Rik; 6030408G03Rik; D130059P03Rik; UBN2

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates using UBN2 Rabbit pAb (A10516) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.

