

KLF15 Rabbit pAb

Catalog NO.: A13420 **1 Publications**

Basic Information

Observed MW

Calculated MW

44kDa

Category

Polyclonal Antibodies

Applications

IF/ICC, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Enables DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in negative regulation of peptidyl-lysine acetylation and positive regulation of transcription by RNA polymerase II. Located in nuclear speck.

Recommended Dilutions

IF/ICC 1:50 - 1:200

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
28999	Q9UIH9

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

KKLF; KLF15

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com


Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Immunofluorescence analysis of A549 cells using KLF15 Rabbit pAb (A13420). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.

