

RNF152 Rabbit pAb

Catalog NO.: A16608 **1 Publications**

Basic Information

Observed MW

22kDa

Calculated MW

22kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Enables small GTPase binding activity and ubiquitin protein ligase activity. Involved in several processes, including cellular response to amino acid starvation; negative regulation of TORC1 signaling; and protein polyubiquitination. Located in lysosomal membrane. Is integral component of organelle membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
220441	Q8N8N0

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

RNF152

Contact

Product Information

☐		order@abclonal.com
☐		service@abclonal.com
☐		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates using RNF152 Rabbit pAb (A16608) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.

