

Pan Lactic acid-Lysine Rabbit pAb

Catalog NO.: A18831SP **1 Publications**

Basic Information

Observed MW

above 10 kDa

Calculated MW

Category

Other Modified Antibodies

Applications

WB, IP, ELISA

Cross-Reactivity

Human, Mouse, Rat, Other (Wide Range)

Conjugation

Unconjugated

Background

Histone lysine lactation (Kla) is a newly discovered histone modification that regulates gene expression in macrophages. In M1 macrophages, lactic acid is derived from incompletely oxidized glucose and then produces lactyl-CoA, which is transferred via acetyltransferase p300 to the lysine tail of the histone. This modification is abundant in gene promoter regions that lack acetylation and are associated with gene expression activation.

Recommended Dilutions

WB 1:500 - 1:1000

IP

0.5 µg - 4 µg antibody for 200 µg - 600 µg extracts of whole cells

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. For high-ratio antibody dilutions (≥1:10000)

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

Product Information

		order@abclonal.com
		service@abclonal.com
		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS, pH 7.3, containing 50% glycerol. Preserved with Proclin300 or sodium azide. May contain 0.05% BSA as specified on the Certificate of Analysis.

Validation Data



Western blot analysis of various lysates using Pan Lactic acid-Lysine Rabbit pAb (A18831SP) at 1:1000 dilution incubated overnight at 4°C. MCF7 cells were treated with Sodium lactate(200 mM) for 24 hours, NIH/3T3 cells were treated with Sodium lactate(100 mM) for 24 hours.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 30 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90 s.

