

OsMADS15 Rabbit pAb

Catalog NO.: A20662

Basic Information

Observed MW

30kDa

Calculated MW

30kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Oryza sativa

Conjugation

Unconjugated

Background

an ortholog of Arabidopsis AP1 gene

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

4342214

Swiss Prot

Q6Q9I2

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

RAP1A; MADS15; FDRMADS3; OsMADS15; RMADS215; OJ1567_G09.109-1

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using OsMADS15 Rabbit pAb (A20662) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 60s.

