

ABflo® 594 Rabbit anti-Mouse/Rat CD31/PECAM1 mAb

Catalog NO.: A23703

Basic Information

Observed MW

Calculated MW

69kDa/77kDa/80kDa/81kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse, Rat

Conjugation

ABflo® 594. Ex:588nm. Em:604nm.

Clone number

ARC61018

Background

Predicted to enable protein homodimerization activity; protein phosphatase binding activity; and transmembrane signaling receptor activity. Acts upstream of or within several processes, including Rho protein signal transduction; endothelial cell morphogenesis; and positive regulation of tyrosine phosphorylation of STAT protein. Located in several cellular components, including cell-cell contact zone; external side of plasma membrane; and smooth muscle contractile fiber. Is expressed in several structures, including alimentary system; cardiovascular system; central nervous system; extraembryonic component; and genitourinary system. Human ortholog(s) of this gene implicated in coronary artery disease (multiple); lung non-small cell carcinoma; neuroblastoma; and psoriatic arthritis. Orthologous to human PECAM1 (platelet and endothelial cell adhesion molecule 1).

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

18613

Swiss Prot

Q08481

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Cd31; Pecam; PECAM-1


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.		

Validation Data



Flow cytometry:1X10^6 C57BL/6 mouse splenocytes were surface-stained with ABflo® 594 Rabbit IgG isotype control(A23821,5 µl/Test,left) or ABflo® 594 Rabbit anti-Mouse CD31 mAb(A23703,5 µl/Test,right).



Flow cytometry: 1X10^6 Rat splenocytes were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821,5 µl/Test,left) or ABflo® 594 Rabbit anti-Mouse CD31 mAb (A23703,5 µl/Test,right).

