

APC Rabbit anti-Human CD79a mAb

Catalog NO.: A24200

Basic Information

Observed MW

Refer to figures

Calculated MW

20kDa/25kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Human

Conjugation

APC. Ex:650nm. Em:660nm.

Clone number

ARC51136

Background

The B lymphocyte antigen receptor is a multimeric complex that includes the antigen-specific component, surface immunoglobulin (Ig). Surface Ig non-covalently associates with two other proteins, Ig-alpha and Ig-beta, which are necessary for expression and function of the B-cell antigen receptor. This gene encodes the Ig-alpha protein of the B-cell antigen component. Alternatively spliced transcript variants encoding different isoforms have been described.

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

973

Swiss Prot

P11912

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

IGA; MB1; MB-1; IGAalpha


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.		

Validation Data



Flow cytometry:1X10⁶ Human PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD19 mAb (A27286,5 µl/Test),and then intracellularly-stained with APC Rabbit IgG isotype control (A24173,5µl/Test,left) or APC Rabbit anti-Human CD79a mAb (A24200,5 µl/Test,right). Cells in the lymphocyte gate were used for analysis.

