

ABflo® 647 Rabbit anti-Human CD26/DPP4 mAb

Catalog NO.: A24411

Basic Information

Observed MW

Refer to figures

Calculated MW

88kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

ABflo® 647. Ex:648nm. Em:664nm.

Clone number

ARC59223

Background

The protein encoded by this gene is identical to adenosine deaminase complexing protein-2, and to the T-cell activation antigen CD26. It is an intrinsic membrane glycoprotein and a serine exopeptidase that cleaves X-proline dipeptides from the N-terminus of polypeptides.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

1803

Swiss Prot

P27487

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

DPP4; ADABP; ADCP2; CD26; DPPIV; TP103; dipeptidyl peptidase 4


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.		

Validation Data



Flow cytometry:1X10^6 K-562 cells (negative control,left) and NCI-H2452 cells (right) were surface-stained with ABflo® 647 Rabbit anti-Human CD26/DPP4 mAb(A24411,5 µl/Test,orange line) or ABflo® 647 Rabbit IgG isotype control (A22070,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).

