

ABflo® 594 Rabbit anti-Mouse/Rat CD54/ICAM-1 mAb

Catalog NO.: A24896

Basic Information

Observed MW

Calculated MW

59kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse, Rat

Conjugation

ABflo® 594. Ex:588nm. Em:604nm.

Clone number

ARC63603

Background

This gene encodes an integral membrane protein that binds leukocyte adhesion protein LFA-1. It participates in the innate immune response.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

15894

Swiss Prot

P13597

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CD54; Ly-47; Icam-1; MALA-2

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.


Validation Data



Flow cytometry: 1×10^6 C2C12 cells (negative control, left) and Mouse thymocytes (right) were surface-stained with ABflo® 594 Rabbit anti-Mouse/Rat CD54 mAb (A24896, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).




Flow cytometry: 1×10^6 Rat splenocytes were surface-stained with ABflo® 594 Rabbit anti-Mouse/Rat CD54 mAb (A24896, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained rat spleen cells were used as blank control (red line).



Flow cytometry: 1×10^6 Mouse thymocytes were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Mouse/Rat CD54 mAb (A24896, 5 μ l/Test, right).

Validation Data



Flow cytometry: 1X10⁶ Rat splenocytes were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821,5 µl/Test,left) or ABflo® 594 Rabbit anti-Mouse/Rat CD54 mAb (A24896,5 µl/Test,right).

