

ABflo® 647 Mouse IgG 2a isotype control

Catalog NO.: A25262

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Conjugation

ABflo® 647. Ex:648nm. Em:664nm.

Clone number

AMC50006

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID Swiss Prot

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,pH7.3.

Validation Data



Flow cytometry: 1X10⁶ Raji cells (negative control, left) and Human PBMC (right) were surface-stained with ABflo® 647 Mouse anti-Human CD3 mAb (A25188, 5 µl/Test, orange line) or ABflo® 647 Mouse IgG 2a isotype control (A25262, 5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

