

ABflo® 594 Mouse IgG 1 isotype control

Catalog NO.: A25489

Basic Information

Observed MW**Calculated MW**

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Conjugation

ABflo® 594. Ex:588nm. Em:604nm.

Clone number

AMC50005

Background

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Mouse	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.03% proclin300,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ MCF7 cells (negative control,left) and SH-SY5Y cells (right) were surface-stained with ABflo® 594 Mouse Anti-Human CD56 mAb (5 µl/Test,orange line) or ABflo® 594 Mouse IgG 1 isotype control (A25489,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ SH-SY5Y cells were surface-stained with ABflo® 594 Mouse IgG 1 isotype control (A25489,5 µl/Test,left) or ABflo® 594 Mouse Anti-Human CD56 mAb (5 µl/Test,right).

