

ABflo® 647 Mouse IgG 2b isotype control

Catalog NO.: A25491

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 647. Ex:648nm. Em:664nm.

Clone number

AMC50008

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Mouse	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.03% proclin300,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ HuT 78 cells (negative control,left) and Raji cells (right) were surface-stained with ABflo® 647 Mouse Anti-Human CD45RA mAb (5 µl/Test,orange line) or ABflo® 647 Mouse IgG 2b isotype control (A25491,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ Raji cells were surface-stained with ABflo® 647 Mouse IgG 2b isotype control (A25491,5 µl/Test,left) or ABflo® 647 Mouse Anti-Human CD45RA mAb (5 µl/Test,right).

