

ABflo®450 anti-Mouse CD11c mAb

Catalog NO.: A25657

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 450. Ex:406nm. Em:445nm.

Clone number

ACC0003

Background

Predicted to enable receptor tyrosine kinase binding activity. Involved in positive regulation of gene expression and positive regulation of myelination. Acts upstream of or within defense response to virus. Located in external side of plasma membrane. Is expressed in heart and hemolymphoid system. Orthologous to human ITGAX (integrin subunit alpha X).

Recommended Dilutions

FC

The suggested use of this reagent is $\leq 0.25 \mu\text{g}$ per 10^6 cells in 100 μl volume.

Immunogen Information

Gene ID

16411

Swiss Prot

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Cr4; N418; Cd11c

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Armenian Hamster

Isotype

Armenian hamster / IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer:Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA

