

ABflo® 488 Rabbit anti-Mouse CD107b/LAMP2 mAb

Catalog NO.: A25930

Basic Information

Observed MW

Calculated MW

46kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 488. Ex:491nm. Em:516nm.

Clone number

ARC63617

Background

Enables protein domain specific binding activity. Involved in several processes, including autophagosome maturation; protein stabilization; and protein targeting to lysosome involved in chaperone-mediated autophagy. Acts upstream of or within muscle cell cellular homeostasis. Located in late endosome membrane; lysosomal membrane; and phagocytic vesicle membrane. Is integral component of autophagosome membrane. Is expressed in central nervous system; egg cylinder; lung; and retina. Used to study Danon disease. Human ortholog(s) of this gene implicated in Danon disease and hypertrophic cardiomyopathy. Orthologous to human LAMP2 (lysosomal associated membrane protein 2).

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

16784

Swiss Prot

P17047

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Mac3; LGP-B; CD107b; Lamp-2; Lamp II; Lamp-2a; Lamp-2b; Lamp-2c

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Flow cytometry: 1X10⁶ NIH/3T3-CD107b-shRNA4 cells (negative control, left) and NIH/3T3-scramble cells (right) were surface-stained with ABflo® 488 Rabbit anti-Mouse CD107b mAb (A25930, 5 µl/Test, orange line) or ABflo® 488 Rabbit IgG isotype control (A22069, 5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

