

PE Rabbit anti-Mouse FOXP3 mAb

Catalog NO.: A26227 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

47kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC68296

Background

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene result in the scurfy phenotype (sf). Alternative splicing results in multiple transcript variants.

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

20371

Swiss Prot

Q99JB6

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

sf; JM2; scurfin


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry:1X10⁶ C57BL/6 mouse splenocytes were surface-stained with ABflo® 647 Rabbit anti-Mouse CD4 mAb (A24721,5 µl/Test) and then intracellularly-stained with PE Rabbit IgG isotype control (A24172,5 µl/Test,left) or PE Rabbit anti-Mouse FOXP3 mAb (A26227,5 µl/Test,right).

