

ABflo® 488 Rabbit anti-Mouse Ki67 mAb

Catalog NO.: A26239 **2 Publications**

Basic Information

Observed MW

Calculated MW

351kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

ABflo® 488. Ex:491nm. Em:516nm.

Clone number

ARC69121

Background

Predicted to enable protein C-terminus binding activity. Involved in regulation of chromatin organization. Acts upstream of or within cholangiocyte proliferation; hepatocyte proliferation; and meiotic cell cycle. Located in chromosome; cytoplasm; and nucl

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

17345

Swiss Prot

E9PVX6

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Ki67; Ki-67; D630048A14Rik


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse splenocytes (untreated, left) and C57BL/6 mouse splenocytes (treated with 2.5 µg/mL Con A for 3 days, right) were surface-stained with ABflo® 700 Rabbit anti-Mouse CD19 mAb (A26193,5 µl/Test) and then intracellularly-stained with ABflo® 488 Rabbit anti-Mouse Ki67 mAb (A26239,5 µl/Test).

