

ABflo® 647 Rabbit anti-Mouse Ki67 mAb

Catalog NO.: A26240 **1 Publications**

Basic Information

Observed MW

Calculated MW

351kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

ABflo® 647. Ex:648nm. Em:664nm.

Clone number

ARC69121

Background

Predicted to enable protein C-terminus binding activity. Involved in regulation of chromatin organization. Acts upstream of or within cholangiocyte proliferation; hepatocyte proliferation; and meiotic cell cycle. Located in chromosome; cytoplasm; and nucl

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

17345

Swiss Prot

E9PVX6

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Ki67; Ki-67; D630048A14Rik


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse splenocytes (Low Expression,left) and RAW 264.7 cells (right) were intracellularly-stained with ABflo® 647 Rabbit anti-Mouse Ki67 mAb (A26240,5 µl/Test,orange line) or ABflo® 647 Rabbit IgG isotype control (A22070,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ RAW 264.7 cells were intracellularly-stained with ABflo® 647 Rabbit IgG isotype control (A22070,5 µl/Test,left) or ABflo® 647 Rabbit anti-Mouse Ki67 mAb (A26240,5 µl/Test,right).

