

ABflo® 488 Mouse IgG 2c isotype control

Catalog NO.: A26279 **1 Publications**

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Conjugation

ABflo® 488. Ex:491nm. Em:516nm.

Clone number

AMC50032_M1

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Mouse	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ 293T cells (negative control,left) and 293T (Transfection,right) cells were surface-stained with Mouse anti-Human CXCR2/CD182 mAb (5 µl/Test,orange line) , followed by Alexa Fluor® 488 conjugated goat anti-mouse pAb staining, or ABflo® 488 Mouse IgG 2c isotype control(A26279,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).

