

APC Rabbit anti-Human IgG (Fc) mAb

Catalog NO.: A26395 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

36kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

APC. Ex:650nm. Em:660nm.

Clone number

ARC53383

Background

Predicted to enable antigen binding activity and immunoglobulin receptor binding activity. Predicted to be involved in several processes, including activation of immune response; defense response to other organism; and phagocytosis. Predicted to act upstream of or within several processes, including immunoglobulin mediated immune response; positive regulation of hypersensitivity; and positive regulation of phagocytosis. Located in extracellular space.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

3500

Swiss Prot

P01857

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

IGHG1


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry:1X10⁶ Human PBMC were intracellularly-stained with FITC Mouse anti-Human CD19 mAb (A22815,5 µl/Test) and APC Rabbit IgG isotype control (A24173,5 µl/Test,left) or APC Rabbit anti-Human IgG (Fc) mAb (A26395,5 µl/Test,right).

