

PE Rabbit anti-Human MRC2 mAb

Catalog NO.: A26561 **1 Publications**

Basic Information

Observed MW

Calculated MW

167kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC66230

Background

This gene encodes a member of the mannose receptor family of proteins that contain a fibronectin type II domain and multiple C-type lectin-like domains. The encoded protein plays a role in extracellular matrix remodeling by mediating the internalization and lysosomal degradation of collagen ligands. Expression of this gene may play a role in the tumorigenesis and metastasis of several malignancies including breast cancer, gliomas and metastatic bone disease.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

9902

Swiss Prot

Q9UBG0

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CD280; UPARAP; CLEC13E; ENDO180


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ Daudi cells (negative control,left) and U-138MG cells (right) were surface-stained with PE Rabbit anti-Human MRC2 mAb (A26561,5 µl/Test,orange line) or PE Rabbit IgG isotype control (A24172,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ U-138MG cells were surface-stained with PE Rabbit IgG isotype control (A24172,5 µl/Test,left) or PE Rabbit anti-Human MRC2 mAb (A26561,5 µl/Test,right).

