

ABflo® 610 Rabbit anti-Human CD107b/LAMP2 mAb

Catalog NO.: A26571 **1 Publications**

Basic Information

Observed MW

Calculated MW

44kDa/45kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Human

Conjugation

ABflo® 610. Ex:421nm. Em:612nm.

Clone number

ARC54767

Background

The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins.

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

3920

Swiss Prot

P13473

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

DND; LAMPB; CD107b; LAMP-2; LGP-96; LGP110


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ knockout (KO) HeLa cells (negative control, left) and HeLa cells (right) were intracellularly-stained with ABflo[®] 610 Rabbit anti-Human LAMP2 mAb (A26571, 5 µl/Test, orange line) or ABflo[®] 610 Rabbit IgG isotype control (A25826, 5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

