

PE Rabbit anti-Human CD62E/E-Selectin mAb

Catalog NO.: A26630 **1 Publications**

Basic Information

Observed MW

Calculated MW

66kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC61987

Background

The protein encoded by this gene is found in cytokine-stimulated endothelial cells and is thought to be responsible for the accumulation of blood leukocytes at sites of inflammation by mediating the adhesion of cells to the vascular lining. It exhibits structural features such as the presence of lectin- and EGF-like domains followed by short consensus repeat (SCR) domains that contain 6 conserved cysteine residues. These proteins are part of the selectin family of cell adhesion molecules. Adhesion molecules participate in the interaction between leukocytes and the endothelium and appear to be involved in the pathogenesis of atherosclerosis.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

6401

Swiss Prot

P16581

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

ELAM; ESEL; CD62E; ELAM1; LECAM2; selectin-e

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1X10⁶ HUVEC cells (untreated,left) and HUVEC cells (treated with 10ng/ml TNF-α for 6 hours,right) were surface-stained with PE Rabbit anti-Human CD62E/E-Selectin mAb (A26630,5 μl/Test,orange line) or PE Rabbit IgG isotype control (A24172,5 μl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ HUVEC cells (untreated,left) and HUVEC cells (treated with 10ng/ml TNF-α for 6 hours,right) were surface-stained with PE Rabbit anti-Human CD62E/E-Selectin mAb (A26630,5 μl/Test).

