

PE Rabbit anti-Mouse Ly-49D mAb

Catalog NO.: A26717 **1 Publications**

Basic Information

Observed MW

Calculated MW

31kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC69648

Background

Predicted to enable carbohydrate binding activity. Predicted to act upstream of or within cell adhesion. Located in external side of plasma membrane.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

16635

Swiss Prot

Q60651

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Chok; Ly49d; Klra32; Klra33; Ly-49d; ly49r<129>


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ 293T cells (negative control,left) and 293T (Transfection,right) cells were surface-stained with PE Rabbit anti-Mouse Ly-49D mAb (A26717,5 µl/Test,orange line) or PE Rabbit IgG isotype control (A24172,5 µl/Test,blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ C57BL/6 mouse splenocytes were surface-stained with PE Rabbit IgG isotype control (A24172,5 µl/Test,left) or PE Rabbit anti-Mouse Ly-49D mAb (A26717,5 µl/Test,right).

