

APC Rat anti-Mouse CD206/MMR mAb

Catalog NO.: A26785 **3 Publications**

Basic Information

Observed MW

Calculated MW

165kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

APC. Ex:650nm. Em:660nm.

Clone number

ATC0035

Background

Enables mannose binding activity and transmembrane signaling receptor activity. Acts upstream of or within cellular response to interferon-gamma; cellular response to interleukin-4; and cellular response to lipopolysaccharide.

Located in cell surface. Is expressed in several structures, including hemolymphoid system; intestine; liver; lung; and yolk sac. Orthologous to human MRC1 (mannose receptor C-type 1).

Recommended Dilutions

FC (intra)

The suggested use of this reagent is $\leq 0.25 \mu\text{g}$ per 10^6 cells in 100 μl volume.

Immunogen Information

Gene ID

17533

Swiss Prot

Q61830

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

MR; CD206

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rat

Isotype

Rat / IgG2b, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.09% sodium azide,pH7.2.

