

APC Rabbit anti-Mouse KLRG1/MAFA mAb

Catalog NO.: A27186

Basic Information

Observed MW

Calculated MW

21kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

APC. Ex:650nm. Em:660nm.

Clone number

ARC67806

Background

Predicted to enable carbohydrate binding activity. Predicted to act upstream of or within cell surface receptor signaling pathway. Predicted to be located in intracellular membrane-bounded organelle. Orthologous to human KLRG1 (killer cell lectin like receptor G1).

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

50928

Swiss Prot

O88713

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

MAFA; 2F1-Ag; MAFA-L


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry:1X10⁶ C57BL/6 mouse splenocytes were surface-stained with ABflo® 488 Rabbit anti-Mouse CD161c/NK1.1 mAb (A24922,5 µl/Test) and APC Rabbit IgG isotype control (A24173,5 µl/Test,left) or APC Rabbit anti-Mouse KLRG1/MAFA mAb (A27186,5 µl/Test,right). Cells in the lymphocyte gate were used for analysis.

