

PE Rabbit anti-Human CD336/NKp44 mAb

Catalog NO.: A27276

Basic Information

Observed MW

Calculated MW

31kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC69299

Background

Predicted to enable signaling receptor activity. Predicted to be involved in cellular defense response and signal transduction. Predicted to be located in plasma membrane. Predicted to be integral component of plasma membrane. Predicted to be active in cell surface.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

9436

Swiss Prot

O95944

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

LY95; CD336; NKP44; NK-p44; dJ149M18.1


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ Human PBMC (treated with rhIL-2 for 7 days) were surface-stained with APC Mouse anti-Human CD56 mAb (A26869, 5 µl/Test) and PE Rabbit IgG isotype control (A24172,5 µl/Test,left) or PE Rabbit anti-Human CD336/NK-p44 mAb (A27276,5 µl/Test,right). Cells in the lymphocyte gate were used for analysis.

