

ABflo® 450 Rat anti-Mouse Ly-6C mAb

Catalog NO.: A27623

Basic Information

Observed MW

Calculated MW

14KDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 450. Ex:406nm. Em:445nm.

Clone number

ATC0057

Background

Predicted to enable acetylcholine receptor binding activity and acetylcholine receptor inhibitor activity. Predicted to be involved in acetylcholine receptor signaling pathway. Located in external side of plasma membrane. Is expressed in embryo mesenchyme; gut; liver; metanephros; and nasal cavity. Orthologous to human LY6H (lymphocyte antigen 6 family member H).

Recommended Dilutions

FC ≤0.25 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

17067

Swiss Prot

P0CW02

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Ly6c; Ly-6C; Ly-6C1

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rat

Isotype

Rat / IgG2c, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,pH7.2.

