

Magnetic Beads-conjugated Rabbit anti-Human TMEM192 mAb [MJFF-G-20C1]

Catalog NO.: A27784

Basic Information

Observed MW

31 kDa

Calculated MW

31 kDa

Category

Monoclonal Antibodies

Applications

IP

Cross-Reactivity

Human

Conjugation

Magnetic Beads

Clone number

ARC5176-Magnetic Beads

Background

Enables protein homodimerization activity. Located in several cellular components, including late endosome; lysosomal membrane; and perinuclear region of cytoplasm.

Recommended Dilutions

IP

20 µL antibody (bead slurry) for 200 µg-600µg extracts of whole cells

Immunogen Information

Gene ID

201931

Swiss Prot

Q8IY95

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

TMEM192


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze / thaw cycles.		
Buffer: PBS with 0.09% Sodium azide,pH7.3.		

Validation Data



Immunoprecipitation of Human TMEM192 from 600 µg extracts of U-251 MG cells was performed using 20 uL of Magnetic Beads-conjugated Rabbit anti-Human TMEM192 mAb [MJFF-G-20C1](A27784). Magnetic Beads-conjugated Rabbit IgG isotype control (AC047) was used to precipitate the Control IgG sample. IP samples were eluted with 1x Laemmli Buffer. The Input lane represents 10% of the total input. Western blot analysis of immunoprecipitates was conducted using TMEM192 Rabbit mAb (A21165) at a dilution of 1:1000.

