

ABflo® 450 Rabbit anti-Human/Monkey CD64/FcγRI mAb

Catalog NO.: A27886

Basic Information

Observed MW

Calculated MW

43kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human, Cynomolgus

Conjugation

ABflo® 450. Ex:406nm. Em:445nm.

Clone number

ARC53680

Background

This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

Hu 2209 /Cyno 102147198

Swiss Prot

P12314

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CD64; FCG1; FCRI; CD64A; FCGR1; IGFR1; FcgammaRI


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with ABflo® 500 Rabbit anti-Human/Monkey CD14 mAb (A26871, 5 µl/Test) and ABflo® 450 Rabbit IgG isotype control (A27451, 5 µl/Test, left) or ABflo® 450 Rabbit anti-Human/Monkey CD64 mAb (A27886, 5 µl/Test, right). Cells in the lymphocyte and monocyte gates were used for analysis.



Flow cytometry: 1X10⁶ Cynomolgus peripheral blood mononuclear cells were surface-stained with ABflo® 450 Rabbit IgG isotype control (A27451, 5 µl/Test, left) or ABflo® 450 Rabbit anti-Human/Monkey CD64/FcγRI mAb (A27886, 5 µl/Test, right). Cells in the monocyte gate were used for analysis.

