

APC/Cyanine7 Rabbit anti-NHP IgD mAb

Catalog NO.: A27949

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Cynomolgus, Rhesus

Conjugation

APC-Cy7. Ex:651nm. Em:779nm.

Clone number

ARC60233

Background

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		



Validation Data



Flow cytometry: 1X10⁶ Cynomolgus PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD19 mAb (A27286,5 µl/Test) and APC/Cyanine7 Rabbit IgG isotype control (5 µl/Test,left) or APC/Cyanine7 Rabbit anti-NHP IgD mAb (A27949,5 µl/Test,right). Cells in the lymphocyte gate were used for analysis.



Flow cytometry: 1X10⁶ Rhesus PBMC were surface-stained with ABflo® 450 Mouse anti-Human CD11b mAb (A27475,5 µl/Test), ABflo® 488 Rabbit anti-Human/Monkey CD19 mAb (A23008,5 µl/Test) and APC/Cyanine7 Rabbit IgG isotype control (5 µl/Test,left) or APC/Cyanine7 Rabbit anti-NHP IgD mAb (A27949,5 µl/Test,right). Cells in the CD11b- gate were used for analysis.



Flow cytometry: 1×10^6 Human peripheral blood mononuclear cells (negative control, left) and Cynomolgus peripheral blood mononuclear cells (right) were surface-stained with APC/Cyanine7 Rabbit anti-NHP IgD mAb (A27949, 5 μ l/Test, right). Cells in the lymphocyte gate were used for analysis.