

ABflo® 610 Rabbit anti-Mouse Ly-6G mAb

Catalog NO.: A28201

Basic Information

Observed MW

Calculated MW

28kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 610. Ex:421nm. Em:612nm.

Clone number

ARC75413

Background

Located in external side of plasma membrane.

Recommended Dilutions

FC ≤0.25 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

546644

Swiss Prot

P35461

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Gr1; Gr-1; Ly-6G


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10^6 C57BL/6 mouse bone marrow cells were surface-stained with ABflo® 450 Rabbit anti-Mouse CD11b mAb (A26229,5 µl/Test) and ABflo® 610 Rabbit IgG isotype control (5 µl/Test,left) or ABflo® 610 Rabbit anti-Mouse Ly-6G mAb (A28201,0.25 µg,right).

