

PerCP Rabbit anti-Human/Mouse FOXP3 mAb

Catalog NO.: A28271

Basic Information

Observed MW

Refer to figures

Calculated MW

47 kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Human, Mouse

Conjugation

PerCP. Ex:482nm. Em:678nm.

Clone number

ARC74142

Background

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

FC (intra) 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

Hu 50943 Me 20371 Cyon 102120882

Swiss Prot

Q9BZS1

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

JM2; AIID; IPEX; PIDX; XPID; DIETER

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1×10^6 Human PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD4 mAb (A26523, 5 μ l/Test) , and then intracellularly-stained with PerCP Rabbit IgG isotype control (A24204, 5 μ l/Test, left) or PerCP Rabbit anti-Human/Mouse/Monkey FOXP3 mAb (A28271, 5 μ l/Test, right). Cells in the lymphocyte gate were used for analysis.



Flow cytometry: 1×10^6 C57BL/6 mouse splenocytes were surface-stained with ABflo® 450 Rabbit anti-Mouse CD4 mAb (A25941, 5 μ l/Test) , and then intracellularly-stained with PerCP Rabbit IgG isotype control (A24204, 5 μ l/Test, left) or PerCP Rabbit anti-Human/Mouse/Monkey FOXP3 mAb (A28271, 5 μ l/Test, right). Cells in the lymphocyte gate were used for analysis.

