

APC Mouse anti-Human CD300e mAb

Catalog NO.: A28307

Basic Information

Observed MW

Calculated MW

23kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

APC. Ex:650nm. Em:660nm.

Background

This gene encodes a member of the CD300 glycoprotein family of cell surface proteins expressed on myeloid cells. The protein interacts with the TYRO protein tyrosine kinase-binding protein and is thought to act as an activating receptor.

Recommended Dilutions

FC ≤0.25 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

342510

Swiss Prot

Q496F6

Immunogen

This information is considered to be commercially sensitive.

Synonyms

CLM2; CLM-2; IREM2; PiGR2; IREM-2; PiGR-2; CD300LE; CMRF35-A5

Contact

Product Information

		order@abclonal.com
		service@abclonal.com
		www.abclonal.com

Source

Mouse

Isotype

Mouse / IgG1, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide and BSA,pH7.2.

