

ABflo® 450 Mouse anti-Mouse TCF-7/TCF-1 mAb

Catalog NO.: A28352

Basic Information

Observed MW

Refer to figures

Calculated MW

45 kDa

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

ABflo® 450. Ex:406nm. Em:445nm.

Background

This gene encodes a transcription factor which is a member of the T-cell specific transcription factor family. The encoded protein is distinct from the hepatic transcription factor, transcription factor 1, which is also referred to by the symbol Tcf1. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined.

Recommended Dilutions

FC (intra)

≤0.125 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

21414

Swiss Prot

Q00417

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Tcf1; TCF-1

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Mouse

Isotype

Mouse IgG1, κ

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide and BSA,pH7.2.

