

APC/Cyanine7 Rat anti-Mouse CD62L/L-Selectin mAb

Catalog NO.: A28368

Basic Information

Observed MW

Calculated MW

42 kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

APC-Cy7. Ex:651nm. Em:779nm.

Clone number

ATC10014

Background

Enables cell adhesion molecule binding activity. Acts upstream of or within response to ATP. Located in external side of plasma membrane. Is expressed in brain; brainstem; submandibular gland epithelium; and testis. Used to study type 1 diabetes mellitus. Human ortholog(s) of this gene implicated in Crohn's disease; IgA glomerulonephritis; and ulcerative colitis. Orthologous to human SELL (selectin L).

Recommended Dilutions

FC ≤0.125 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

20343

Swiss Prot

P18337

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Lnhr; CD62L; LAM-1; Ly-22; Lyam1; LECAM1; Ly-m22; Lyam-1; LECAM-1; L-selectin

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rat

Isotype

Rat / IgG2a, kappa

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,pH7.2.

