

Rabbit anti-Human CD52 mAb

Catalog NO.: A28616

Basic Information

Observed MW

Calculated MW

7kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

Unconjugated

Clone number

ARC69477

Background

Involved in positive regulation of cytosolic calcium ion concentration. Predicted to be located in extracellular region and plasma membrane. Predicted to be intrinsic component of plasma membrane. Predicted to be active in sperm midpiece.

Recommended Dilutions

FC 1:100-1:500

Immunogen Information

Gene ID

1043

Swiss Prot

P31358

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HE5; CDW52; EDDM5


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS,pH7.3.		

Validation Data



Flow cytometry: 1X10^6 Human PBMC were surface-stained with Rabbit anti-Human CD52 mAb (A28616,2 µg/mL,orange line) or Rabbit IgG isotype control (AC042,2 µg/mL,blue line), followed by Alexa Fluor® 647 conjugated goat anti-rabbit pAb staining. Non-fluorescently stained cells were used as blank control (red line). Cells in the lymphocyte gate were used for analysis.

