

PerCP/Cyanine5.5 Rabbit anti-NHP IgD mAb

Catalog NO.: A28680

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Cynomolgus, Rhesus

Conjugation

PerCP-Cy5.5. Ex:482nm. Em:695nm.

Clone number

ARC60233

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

无

Swiss Prot

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit		Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ Human PBMC (negative control,left) and Rhesus PBMC (right) were surface-stained with ABflo[®] 450 Rabbit anti-Human/Monkey CD19 mAb (A27286,5 µl/Test) and PerCP/Cyanine5.5 Rabbit IgG isotype control (A25620,5 µl/Test) or PerCP/Cyanine5.5 Rabbit anti-NHP IgD mAb (A28680,5 µl/Test).



Flow cytometry: Flow cytometry:1X10⁶ Human PBMC (negative control,left) and Cynomolgus PBMC (right) were surface-stained with ABflo[®] 450 Rabbit anti-Human/Monkey CD19 mAb (A27286,5 µl/Test) and PerCP/Cyanine5.5 Rabbit IgG isotype control (A25620,5 µl/Test) or PerCP/Cyanine5.5 Rabbit anti-NHP IgD mAb (A28680,5 µl/Test).

