

PE Rabbit anti-Monkey CD159a mAb

Catalog NO.: A29251

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Cynomolgus

Conjugation

PE. Ex:565nm. Em:574nm.

Clone number

ARC82681

Background

CD159a (also known as NKG2A, encoded by the KLRC1 gene) is an inhibitory receptor present on the surface of natural killer (NK) cells and some T cells. When it binds to the HLA-E molecule, it transmits inhibitory signals and regulates the activity of immune cells. Certain viruses (such as HIV and CMV) upregulate the expression of HLA-E to bind to CD159a and inhibit the function of NK cells. Blocking the NKG2A/HLA-E pathway can restore antiviral immunity.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

102145514

Swiss Prot

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

nkg2A

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1X10⁶ Cynomolgus PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD3 mAb (A27177, 5 µl/Test) and PE Rabbit anti-Monkey CD159a mAb (A29251, 5 µl/Test). Cells in the lymphocyte gate were used for analysis.

