

CLEC3B/Tetranectin Rabbit pAb

Catalog NO.: A6418

Basic Information

Observed MW

Refer to figures

Calculated MW

21 kDa

Category

Polyclonal Antibodies

Applications

WB, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Background

Enables calcium ion binding activity; heparin binding activity; and kringle domain binding activity. Involved in bone mineralization and cellular response to transforming growth factor beta stimulus. Located in cytoplasm; extracellular space; and granular component. Part of collagen-containing extracellular matrix. Implicated in osteoarthritis.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:10 - 1:100

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7123

Swiss Prot

P05452

Immunogen

This information is considered to be commercially sensitive.

Synonyms

TN; TNA; MCDR4; CLEC3B/Tetranectin

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com


Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS, pH 7.3, containing 50% glycerol. Preserved with Proclin300 or sodium azide. May contain 0.05% BSA as specified on the Certificate of Analysis.

Validation Data



Immunofluorescence analysis of HeLa cells using CLEC3B/Tetranectin Rabbit pAb (A6418). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.

