

TFPT Rabbit pAb

Catalog NO.: A6461

Basic Information

Observed MW

Refer to Figures

Calculated MW

28kDa

Category

Polyclonal Antibodies

Applications

WB, IF/ICC, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Predicted to enable DNA binding activity and protein kinase binding activity. Involved in apoptotic signaling pathway. Located in nucleoplasm. Part of Ino80 complex.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:100

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

29844

Swiss Prot

P0C1Z6

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

FB1; amida; INO80F; TFPT

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com


Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Immunofluorescence analysis of U2OS cells using TFPT Rabbit pAb (A6461).Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution.

