

MEPCE Rabbit pAb

Catalog NO.: A7121

Basic Information

Observed MW

Refer to figures

Calculated MW

74kDa

Category

Polyclonal Antibodies

Applications

IF/ICC, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Enables 7SK snRNA binding activity and RNA 5'-methyltransferase activity. Involved in RNA modification; positive regulation of protein localization to Cajal body; and positive regulation of snRNA transcription by RNA polymerase II. Located in nucleus. Part of 7SK snRNP.

Recommended Dilutions

IF/ICC 1:50 - 1:100

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

56257

Swiss Prot

Q7L2J0

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

BCDIN3; MEPCE


Contact

Product Information

☐		order@abclonal.com
☐		service@abclonal.com
☐		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at -20°C. Avoid freeze / thaw cycles.		
Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.		

Validation Data



Immunofluorescence analysis of A549 cells using MEPCE Rabbit pAb (A7121).Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution.

