

Biotin Conjugated-Rabbit IgG isotype control

Catalog NO.: AC042BC

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

IHC-P

Cross-Reactivity

Conjugation

Biotin

Clone number

ARC5105-03

Background

The protein encoded by this gene is a transcriptional regulator and tumor suppressor, serving as an activator of genes involved in both innate and acquired immune responses.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID	Swiss Prot
	无

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human tonsil tissue using Biotin Conjugated-Rabbit IgG isotype control (AC042BC) at a dilution of 1:50 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.

