

ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb

Catalog NO.: AS155

Basic Information

Observed MW

Calculated MW

36kDa

Category

Secondary Antibodies

Applications

FC

Cross-Reactivity

Mouse

Conjugation

ABflo® 488. Ex:491nm. Em:516nm.

Clone number

ARC68493

Background

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

/

Swiss Prot

P01868

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information


☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ Raji cells (negative control,left) and Human PBMC (right) were surface-stained with Mouse anti Human CD3 mAb (2 µg/mL,orange line) or Mouse IgG isotype control (2 µg/mL,blue line), followed by ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb (AS155,5 µl/Test) staining. Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with Mouse IgG isotype control (2 µg/mL,left) or Mouse anti Human CD3 mAb (2 µg/mL,right). ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb (AS155,5 µl/Test) was used as a secondary antibody.

