

pan-Phospho-Serine Rabbit pAb

Catalog No.: AP0932

6 Publications

Basic Information

Observed MW

60kDa

Calculated MW

Category

Polyclonal Antibodies

Applications

WB, DB, ELISA

Cross Reactivity

Mouse

Background

Recommended Dilutions

WB 1:500 - 1:1000

DB 1:1000 - 1:3000

ELISA Recommended starting concentration is 1 µg /mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

Swiss-Prot

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purity

Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of lysates from NIH/3T3 cells, using pan-Phospho-Serine Rabbit pAb (AP0932) at 1:1000 dilution. NIH/3T3 cells were treated with Calyculin A (100 nM) at 37°C for 30 minutes after serum-starvation overnight.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 180s.



Dot-blot analysis of all sorts of peptides using pan-Phospho-Serine Rabbit pAb (AP0932) at 1:1000 dilution.