

THBS2 Rabbit pAb

Catalog No.: A8561**8 Publications**

Basic Information

Observed MW

210KDa

Calculated MW

130 kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross Reactivity

Human

Background

The protein encoded by this gene belongs to the thrombospondin family. It is a disulfide-linked homotrimeric glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein has been shown to function as a potent inhibitor of tumor growth and angiogenesis. Studies of the mouse counterpart suggest that this protein may modulate the cell surface properties of mesenchymal cells and be involved in cell adhesion and migration.

Recommended Dilutions

WB 1:1000 - 1:2000**ELISA**

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7058

Swiss-Prot

P35442

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TSP2; THBS2

Contact

 | 400-999-6126 | cn.market@abclonal.com.cn | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG


Purity

Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of lysates from SK-N-SH cells using THBS2 Rabbit pAb (A8561) at 1:2000 dilution incubated at room temperature for 1.5 hours.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90 s.