

Recombinant Human IFI30 Protein

Catalog No.: RP02995

Sequence Information

Species **Gene ID**
Human 10437

Swiss Prot **Tags**
P13284 C-His

Synonyms
GILT, IFI-30, IFI30, IP-30, IP30

Product Information

Source **Purification**
HEK293 cells ≥ 95 % as
determined by
SDS-PAGE.

Calculated MW **Observed MW**
24.7 kDa 31-34 kDa

Endotoxin
< 1 EU/μg of the protein by LAL
method.

Formulation
Lyophilized from sterile PBS, pH
7.4

Reconstitution
Centrifuge the vial before opening.
Reconstitute to a concentration of
0.1-0.5 mg/mL in sterile distilled
water. Avoid vortex or vigorously
pipetting the protein. For long
term storage, it is recommended
to add a carrier protein or stabilizer
(e.g. 0.1% BSA, 5% HSA, 10% FBS
or 5% Trehalose), and aliquot the
reconstituted protein solution to
minimize free-thaw cycles.

Background

Basic Information

Description

Recombinant Human IFI30 Protein is produced by HEK293 cells expression system. The target protein is expressed with sequence (Met1-Lys243) of human IFI30 (Accession #NP_006323.2) fused with a 6×His tag at the C-terminus.

Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Storage

Store at -20°C

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.


Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Validation Data



Recombinant Human IFI30 Protein was determined by SDS-PAGE under reducing conditions with Coomassie Blue.